




Part 1 Liter Vessel
Part Number 14537-01
Serial Number 101833 - 101875



| | | Tolerances | | |
|------------|--------------------|---------------------------|-------------------|----------|
| Dimensions | Vessel | Inside Ø a | 102.9 to 105.1 mm | -- |
| | | Height b | 186.9 mm | ± 0.7 mm |
| | | With hemispherical bottom | | |
| Material | Borosilicate glass | | | |

Aesch, 2025-01-15

Tayyib Özüdogru
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate