




Part 1 Liter Vessel inactinic brown  
Part Number 14537-02  
Serial Number 100103 - 100104



			Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	--
		Height b	186.9 mm	± 0.7 mm
		With hemispherical bottom		
Material	Borosilicate glass			

Aesch, 2024-07-31

Tayyib Özüdogru  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

# Test Certificate