




Part 1 Liter Vessel inactinic brown  
Part Number 14537-02  
Serial Number 100109 - 100115



|            |                    | Tolerances                |                   |          |
|------------|--------------------|---------------------------|-------------------|----------|
| Dimensions | Vessel             | Inside Ø a                | 102.9 to 105.1 mm | --       |
|            |                    | Height b                  | 186.9 mm          | ± 0.7 mm |
|            |                    | With hemispherical bottom |                   |          |
| Material   | Borosilicate glass |                           |                   |          |

Aesch, 2024-07-31

Tayyib Özüdogru  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

# Test Certificate